

Product Datasheet

Vaccinia L1R/Protein L1 antibody (orb1290086)

Description	Mouse monoclonal antibody to Vaccinia L1R/Protein L1
Species/Host	Mouse
Reactivity	Virus
Endotoxins	Please contact with the lab for this information.
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Preservatives	PBS buffer, pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 95%
Expiration Date	12 months from date of receipt.