

Product Datasheet

Vaccinia H3L/Ag35 antibody (orb1290087)

Description	Mouse monoclonal antibody to Vaccinia H3L/Ag35
Species/Host	Mouse
Reactivity	Virus
Endotoxins	Please contact with the lab for this information.
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Preservatives	PBS buffer, pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at +4°C short term (1-2 weeks). Store at -20 °C 12 months. Store at -80°C long term.
Note	For research use only
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 95%
Expiration Date	12 months from date of receipt.