

Product Datasheet

Monkeypox virus D14L Protein (orb1290088)

| | |
|----------------------------|--|
| Description | Recombinant Monkeypox virus D14L Protein |
| Reactivity | Virus |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, SDS-PAGE, WB |
| Preservatives | Supplied as solution form in PBS pH 7.4. |
| Form/Appearance | Lyophilized |
| Storage | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt. |
| Note | For research use only |
| Application notes | Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details. |
| Purity | > 90% as determined by SDS-PAGE. |
| MW | 23.73 kDa |
| Expiration Date | 6 months from date of receipt. |