

Product Datasheet

Monkeypox virus B6R/SL-159 Protein (orb1290097)

Description	Recombinant Monkeypox virus B6R/SL-159 Protein
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Preservatives	Supplied as solution form in PBS pH 7.4.
Form/Appearance	Lyophilized
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
Application notes	Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.
Purity	> 90% as determined by SDS-PAGE
MW	31.73 kDa
Expiration Date	6 months from date of receipt.