

## Product Datasheet

# Monkeypox virus B2R/Hemagglutinin Protein (orb1290107)

<b>Description</b>	Recombinant Monkeypox virus B2R/Hemagglutinin Protein
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, SDS-PAGE, WB
<b>Preservatives</b>	Supplied as solution form in PBS pH 7.4.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	Reconstitute in sterile water for a stock solution.A copy of datasheet will be provided with the products, please refer to it for details.
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	30.23 kDa
<b>Expiration Date</b>	6 months from date of receipt.