

Product Datasheet

SARS-CoV-2(2019-nCoV) Spike glycoprotein antibody (orb1290116)

Description	Mouse monoclonal antibody to SARS-CoV-2(2019-nCoV) Spike glycoprotein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 (2019-nCoV) Spike protein fragment 1[679-833]
Preservatives	PBS, pH 7.4, containing 0.05% proclin300, 50% glycerol.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	WB:1:1000-4000
Expiration Date	12 months from date of receipt.