

Product Datasheet

SARS-CoV-2(2019-nCoV) S2 antibody (orb1290139)

| | |
|----------------------------|---|
| Description | Human monoclonal antibody to SARS-CoV-2(2019-nCoV) S2 |
| Species/Host | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Recombinant SARS-CoV-2 S2 protein |
| Preservatives | PBS,pH7.5 |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Human IgG |
| Clonality | Monoclonal |
| Dilution Range | WB:1:1000-2000; ELISA:1:5000-1:10000 |
| Expiration Date | 12 months from date of receipt. |