

Product Datasheet

SARS-CoV-2(2019-nCoV) S2 antibody (orb1290141)

Description	Human monoclonal antibody to SARS-CoV-2(2019-nCoV) S2
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 S2 protein
Preservatives	PBS,pH7.5
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Human IgG
Clonality	Monoclonal
Dilution Range	WB:1:1000-2000; ELISA:1:5000-1:10000
Expiration Date	12 months from date of receipt.