

Product Datasheet

SARS-CoV-2 (2019-nCoV) S RBD (His-Avi,HEK293) Protein (orb1290151)

Description	Recombinant SARS-CoV-2 (2019-nCoV) S RBD (His-Avi,HEK293) Protein
Reactivity	Virus
Endotoxins	Please contact with the lab for this information
Conjugation	Unconjugated
Preservatives	Supplied as solution form in PBS, pH 7.5 or lyophilized from PBS, pH 7.5
Storage	Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt.
Note	For research use only
MW	29.99 kD
Expiration Date	6 months from date of receipt.