

Product Datasheet

Envafolimab biosimilar antibody (orb1290173)

Description	Envafolimab biosimilar antibody
Species/Host	Mouse
Reactivity	Chimeric
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	CD274/PD-L1/B7-H1
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	(VH-CH2-CH3) ₂
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	2102192-68-5
Expiration Date	12 months from date of receipt.