

Product Datasheet

Cinpanemab biosimilar antibody (orb1290175)

Description	Cinpanemab biosimilar antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SNCA
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1-lambda
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	2094516-02-4
Expiration Date	12 months from date of receipt.