

Product Datasheet

Coltuximab biosimilar antibody (orb1290209)

Description	Coltuximab biosimilar antibody
Species/Host	Mouse
Reactivity	Chimeric
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	CD19
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1-kappa
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	1269764-99-9
Expiration Date	12 months from date of receipt.