

## Product Datasheet

### Margetuximab biosimilar antibody (orb1290213)

<b>Description</b>	Margetuximab biosimilar antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Chimeric
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC
<b>Target</b>	ERBB2/EGFR2/CD340
<b>Preservatives</b>	PBS buffer pH 7.5
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1-kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>Source</b>	CHO cells
<b>CAS Number</b>	1350624-75-7
<b>Expiration Date</b>	12 months from date of receipt.