

Product Datasheet

Guselkumab biosimilar antibody (orb1290215)

Description	Guselkumab biosimilar antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	IL23A
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1-lambda
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	1350289-85-8
Expiration Date	12 months from date of receipt.