

Product Datasheet

Bimagrumab biosimilar antibody (orb1290220)

Description	Bimagrumab biosimilar antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	ACVR2A/ACVR2B
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1-lambda
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	1356922-05-8
Expiration Date	12 months from date of receipt.