

## Product Datasheet

### Nivolumab biosimilar antibody (orb1290228)

<b>Description</b>	Nivolumab biosimilar antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, FC
<b>Target</b>	PDCD1/PD1/CD279
<b>Preservatives</b>	PBS buffer pH 7.5
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG4-kappa
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95%
<b>Source</b>	CHO cells
<b>CAS Number</b>	946414-94-4
<b>Expiration Date</b>	12 months from date of receipt.