

Product Datasheet

Palivizumab biosimilar antibody (orb1290232)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Description | Palivizumab biosimilar antibody |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA |
| Target | RSV |
| Preservatives | PBS buffer pH 7.5 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1-kappa |
| Clonality | Monoclonal |
| Purity | > 95% |
| Source | CHO cells |
| CAS Number | 188039-54-5 |
| Expiration Date | 12 months from date of receipt. |