

Product Datasheet

Enoblituzumab biosimilar antibody (orb1290251)

Description	Enoblituzumab biosimilar antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	CD276/B7-H3
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1-kappa
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	1353485-38-7
Expiration Date	12 months from date of receipt.