

Product Datasheet

Talacotuzumab biosimilar antibody (orb1290257)

Description	Talacotuzumab biosimilar antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	NOTCH2/ NOTCH3
Preservatives	PBS buffer pH 7.5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2-kappa
Clonality	Monoclonal
Purity	> 95%
Source	CHO cells
CAS Number	1826831-79-1
Expiration Date	12 months from date of receipt.