

Product Datasheet

Human IL-36 gamma protein (orb1291242)

Description	The Human IL-36 gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-36 gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SMCKPITGTI NDLNQQVWTL QGQNLVAVPR SDSVTPVTVA VITCKYPEAL EQGRGDPIYL GIQNPEMCLY CEKVGEQPTL QLKEQKIMDL YGQPEPVKPF LFYRAKTGRT STLESVAFPD WFIASSKRDQ PIILTSELGK SYNTAFELNI ND (152)
Purity	98%
Source	Yeast
MW	17.0 kDa
Entrez	56300
Expiration Date	6 months from date of receipt.