

Product Datasheet

Canine IL-36RA protein (orb1291244)

Description	The Canine IL-36 Receptor Antagonist yeast-derived recombinant protein is not tagged and is naturall...
Target	IL-36RA
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	VLSGALCFRM KDAALKVLYL HDNQLLAGGL HAGKVIK GEE ISVVPNRSLD AKLSPVILGV QGGSQCLSCG TGQEPTL KLE PVNIMELYLG AKESKSFTFY RRDTGLTSSF ESAAYPGWFL CTVPEADQPL RLTQLSEDAS WDNPITDFYF QQCD (154)
Purity	98%
Source	Yeast
MW	16.9 kDa
Entrez	611869
Expiration Date	6 months from date of receipt.