

Product Datasheet

Guinea Pig IL-6 protein (orb1291249)

Description The Guinea Pig IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....

Target IL-6

Form/Appearance Lyophilized

Storage -20°C

Note For research use only

Protein Sequence AFPTAQVQHD FTADTTDEMT TAEMTTTMPN
KPTTSASQVF QMFMRVYQAV KELKNEMNKH
NVEKAILNNL DLPKLELDG CFFNGYNWET
CQLKITPGLF KFQTYLQSMQ NKLQNESENK
KAANIYAGIK SLSLFMKSKI NNTEQMEFLS
PTPDATLLEK LETQSQTQML LIAEIVLQRL
EEFLQDSLRA IRKADWEGRN (200)

Purity 98%

Source Yeast

MW 23.0 kDa

Entrez [106027696](#)

Expiration Date 6 months from date of receipt.