

## Product Datasheet

### Ferret APRIL (TNFSF13) protein (orb1291250)

<b>Description</b>	The Ferret APRIL yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
<b>Target</b>	APRIL
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	AALMHKHKQK KKHSVLHLVP INITSKEDSD VTEVMWQPAL KRGRGLEAQQ YVVRVWDAGV YLLYSQVLFH DVTFTMGQVW SREGQGRQET LFRCVRSMP5 NPDWAYNSCY SAGVFHLHQG DALSVTIPRG RAKLSLSPHG TFLGLVKL (148)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	16.6 kDa
<b>Entrez</b>	<b>116578386, 101683970</b>
<b>Expiration Date</b>	6 months from date of receipt.