

## Product Datasheet

### CCL3 antibody (Biotin) (orb1291261)

<b>Description</b>	The Black Flying Fox Bat CCL3 polyclonal antibody is labeled and has been qualified for use in ELISA...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Bat
<b>Conjugation</b>	Biotin
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant Bat CCL3
<b>Target</b>	CCL3
<b>Storage</b>	4 degrees C.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.