

Product Datasheet

CXCL6 antibody (Biotin) (orb1291267)

Description	The Ferret CXCL6 polyclonal antibody is labeled and has been qualified for use in ELISA. It is the d...
Species/Host	Rabbit
Reactivity	Ferret
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Ferret CXCL6
Target	CXCL6
Storage	4 degrees C.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.