

Product Datasheet

Mouse CTLA4 ELISA Kit (orb1291269)

Description	Mouse CTLA4 (Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit
Reactivity	Mouse
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids
Uniprot ID	P09793
Solubility (25°C)	2000 pg/mL
Sensitivity	12.4 pg/mL
Expiration Date	Please enquire.

