

## Product Datasheet

### Human OPI ELISA Kit (orb1291278)

<b>Description</b>	Human OPI (Opiorphin) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.16-10 ng/mL
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive Inhibition
<b>Sample Types</b>	tissue homogenates, cell lysates, cell culture supernates and other biological fluids
<b>Uniprot ID</b>	<a href="#">Q99935</a>
<b>Solubility (25°C)</b>	10 ng/mL
<b>Sensitivity</b>	0.049 ng/mL
<b>Expiration Date</b>	Please enquire.

