

Product Datasheet

General PGE1 ELISA Kit (orb1291328)

Description	General PGE1 (Prostaglandin E1) ELISA Kit
Reactivity	All
Range	31.25-2000 pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Competitive Inhibition
Sample Types	Cell culture supernatant, Plasma, Saliva, Serum, Tissue Culture Media, Urine
Solubility (25°C)	2000 pg/mL
Sensitivity	9.55 pg/mL
Expiration Date	Please enquire.

