

## Product Datasheet

### ELISA kit for Porcine TNF-alpha (ALP) (orb1291397)

<b>Description</b>	ELISA kit for Porcine TNF-alpha (ALP)
<b>Reactivity</b>	Porcine
<b>Range</b>	10-1000 pg/ml
<b>Conjugation</b>	ALP
<b>Tested Applications</b>	ELISA
<b>Form/Appearance</b>	without plates
<b>Storage</b>	18 months
<b>Note</b>	For research use only
<b>Expiration Date</b>	Please enquire.