

Product Datasheet

MPXV A29 antibody (orb1291545)

Description	MPXV A29 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Monkeypox Virus (MPXV) A29 Protein
Target	MPXV A29
Preservatives	10mM PBS(pH7.4),0.03% Proclin300
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	12.5 kDa
Uniprot ID	Q77HM6
Expiration Date	12 months from date of receipt.