



Product Datasheet

GPCR RDC1/CXCR-7/ACKR3 Antibody (orb1291687)


Description	GPCR RDC1/CXCR-7/ACKR3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human GPCR RDC1/CXCR-7/ACKR3 (RNYRYELMKAFIFKYSAKTGLTKLIDASRVSE).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P25106
Expiration Date	12 months from date of receipt.



IHC analysis
of GPCR
RDC1/CXCR-
7/ACKR3 u...



IHC analysis
of GPCR
RDC1/CXCR-
7/ACKR3 u...



IHC analysis
of GPCR
RDC1/CXCR-
7/ACKR3 u...