



Product Datasheet

ZNRF3 Antibody (orb1291688)


Description	ZNRF3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, WB
Immunogen	E.coli-derived human ZNRF3 recombinant protein (Position: R54-E913).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	130 kDa
Uniprot ID	Q9ULT6
Expiration Date	12 months from date of receipt.



IF analysis of ZNRF3 using anti-ZNRF3 an...



Western blot analysis of ZNRF3 using ant...



Flow Cytometry analysis of PC-3 cells us...