




Product Datasheet

UBC12/UBE2M Antibody (orb1291705)

Description	UBC12/UBE2M Antibody	 <p>IHC analysis of UBC12/UBE2M using anti-U...</p>
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, FC, ICC, IF, IHC, WB	 <p>IHC analysis of UBC12/UBE2M using anti-U...</p>
Immunogen	E.coli-derived human UBC12/UBE2M recombinant protein (Position: K25-Q154).	
Form/Appearance	Lyophilized	 <p>IHC analysis of UBC12/UBE2M using anti-U...</p>
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.	
Note	For research use only	
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.	
Isotype	Rabbit IgG	
Clonality	Polyclonal	
MW	21 kDa	
Uniprot ID	P61081	
Expiration Date	12 months from date of receipt.	