

Product Datasheet

Rat OX40 / TNFRSF4 / CD134 Protein (orb1291792)

Description	Rat OX40, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 20 - Pro 210 (Acces...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	OX40
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 22.9 kDa. The protein migrates as 43-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Val 20 - Pro 210
Purity	90%
Source	Rat OX40, His Tag (orb1291792) is expressed from human 293 cells (HEK293). It contains AA Val 20 - Pro 210 (Accession # P15725-1).
MW	22.9 kDa
Uniprot ID	P15725
Expiration Date	6 months from date of receipt.