

Product Datasheet

Rabbit Her2 / ErbB2 Protein (orb1291803)

Description	Rabbit Her2, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 23 - Thr 652 (Ac...
Reactivity	Rabbit
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Her2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 71.0 kDa. The protein migrates as 85-95 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_002719389.1
Purity	95%
Source	Rabbit Her2, His Tag (orb1291803) is expressed from human 293 cells (HEK293). It contains AA Thr 23 - Thr 652 (Accession # G1SZL0-1).
MW	71.0 kDa
Uniprot ID	G1SZL0
NCBI	XP_002719389.1
Expiration Date	6 months from date of receipt.