

Product Datasheet

Goat Anti-Rabbit IgG (H+L) (orb1291811)

Description	Goat Anti-Rabbit IgG (H+L)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Immunogen	Rabbit IgG (H+L)
Form/Appearance	Liquid: 0.01M PBS, pH 7.4.
Storage	Upon receiving, store at -20°C at least. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	95%. Affinity purification
Expiration Date	12 months from date of receipt.