

Product Datasheet

Goat Anti-Human IgG (Fab) antibody (orb1292422)

Description	Rabbit Polyclonal antibody against Goat Anti-Human IgG (Fab) conjugated to HRP
Species/Host	Goat
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, WB
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:2000-1:10000 for western blotting with ECL substrates1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.