

Product Datasheet

Rabbit Anti-Chicken IgY (IgG)(H+L) antibody (orb1292424)

Description	Rabbit Polyclonal antibody against Rabbit Anti-Chicken IgY (IgG)(H+L) conjugated to Biotin
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Biotin
Tested Applications	ELISA, FC, IF, WB
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:100 - 1:1000 for enzyme immunohistochemistry on tissue sections 1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.