

## Product Datasheet

# Goat Anti-Monkey IgG(gamma) antibody (orb1292426)

<b>Description</b>	Rabbit Polyclonal antibody against Goat Anti-Monkey IgG(gamma) conjugated to HRP
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Monkey
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	ELISA, WB
<b>Preservatives</b>	3mg/ml rAlbumin ( IgG-Free, Protease-Free)
<b>Form/Appearance</b>	liquid
<b>Storage</b>	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20°C for 24 months, avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Hazard Information</b>	This product is for research use only, not for diagnostic or therapeutic use.
<b>Dilution Range</b>	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates
<b>Expiration Date</b>	12 months from date of receipt.