

Product Datasheet

Mouse Anti-Rabbit IgG Heavy Chain Specific antibody (orb1292435)

Description	Rabbit Polyclonal antibody against Mouse Anti-Rabbit IgG Heavy Chain Specific conjugated to HRP
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, WB
Preservatives	0.5% rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, 0.1% proclin300, pH7.3. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Dilution Range	1:5000-1:20000 for western blotting with ECL substrates 1:5000-1:20000 for ELISA and Western blotting with chromogenic substrates
Expiration Date	12 months from date of receipt.