

Product Datasheet

TAU (Phospho-S404) antibody (orb1292452)

Description	Rabbit Polyclonal antibody to TAU (Phospho-S404)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S404 of human Tau
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60kDa
Uniprot ID	P10636
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.