

Product Datasheet

RB (Phospho-S807/811) antibody (orb1292457)

Description	Rabbit Polyclonal antibody to RB (Phospho-S807/811)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S807 & S811 of human RB
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	120 kDa
Uniprot ID	P06400
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.