

## Product Datasheet

### NFKB1 Phospho-S337 antibody (orb1292472)

<b>Description</b>	Rabbit Polyclonal antibody to NFKB1 Phospho-S337
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A phospho specific peptide corresponding to residues surrounding S337 of human NFKB1
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.06% PC300 and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>MW</b>	120 kd
<b>Uniprot ID</b>	<a href="#">P19838</a>
<b>Dilution Range</b>	WB: 1:500-1:2000
<b>Expiration Date</b>	12 months from date of receipt.