

## Product Datasheet

### RFP antibody (orb1292489)

<b>Description</b>	Mouse Monoclonal antibody to RFP
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Red fluorescent protein
<b>Form/Appearance</b>	liquid
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥95% as determined by SDS-PAGE
<b>Uniprot ID</b>	<b>Q8ISF8</b>
<b>Dilution Range</b>	WB: 1:1000-1:5000
<b>Expiration Date</b>	12 months from date of receipt.