

Product Datasheet

PKC (pan) (zeta) Phospho-T410 antibody (orb1292499)

Description	Rabbit Polyclonal antibody to PKC (pan) (zeta) Phospho-T410
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	A synthetic phosphorylated peptide around T410 of human PRKCZ
Target	Phospho-pan PKC (T410)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	85 kDa
Uniprot ID	P05129 , P05771 , P17252 , Q02156 , P24723 , P41743 , Q05655 , Q05513 , Q04759
Dilution Range	WB: 1:500 - 1:2000
Expiration Date	12 months from date of receipt.