

Product Datasheet

IL-16 antibody (orb1292516)

Description	Rabbit Polyclonal antibody to IL-16
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	interleukin 16(lymphocyte chemoattractant factor)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 12 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	59,66,141 kDa
Uniprot ID	Q14005
Dilution Range	WB: 1:500 - 1:2000
Expiration Date	12 months from date of receipt.