

Product Datasheet

HNRNPA1 antibody (orb1292559)

Description	Rabbit Polyclonal antibody to HNRNPA1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	heterogeneous nuclear ribonucleoprotein A1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	34kDa, 38kDa
Uniprot ID	P09651
Dilution Range	WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.