

Product Datasheet

2019 nCOV N protein antibody (orb1292565)

Description	Rabbit Polyclonal antibody to 2019 nCOV N protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, GICA, WB
Form/Appearance	liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	ELISA: Use at an assay dependent concentration (Suggestions: Use a concentration of 0.1 - 1 µg/ml)
Expiration Date	12 months from date of receipt.