


# Product Datasheet

## GABA ELISA Kit (orb1292583)

<b>Description</b>	GABA (Gamma-aminobutyric acid) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	GABA
<b>Storage</b>	2 - 8 °C for 12 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.



GABA (Gamma-aminobutyric acid) ELISA Kit

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)