


## Product Datasheet

### Rat Bgn ELISA Kit (orb1292623)

<b>Description</b>	Rat Bgn ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	7.813-500pg/ml
<b>Target</b>	Bgn
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P47853</b>
<b>Sensitivity</b>	4.688pg/ml
<b>Expiration Date</b>	Please enquire.



Rat Bgn  
(Biglycan)  
ELISA Kit